

INCUT



**VECTEC®**

SINGLE USE  
LAPAROSCOPIC DEVICES





10.9



REF: CAB 3000  
« Single use monopolar cable »

## INCUT

### HIGH QUALITY SINGLE USE LAPAROSCOPIC SCISSORS

Reference	Working length	Diameter	Blade cut length	Connector
IC5 330	330 mm	5 mm	17 mm	Radial
IK5 330	330 mm	5 mm	17 mm	Axial
CS 470*	470 mm	5 mm	17 mm	Radial

\* box of 5 individually packaged units



Ergonomic handle allowing a precise surgical gesture.

IK5 330 model developed for gastrointestinal surgery.

IC5 330 model dedicated to gynecological surgery.

Double action curved blades facilitating the dissection.

Monopolar coagulation with extra long electrical insulation up to the active portion of the blades.

Packaging: box of 25 individually packaged units.

MD Class IIb  
Sterilized by ETO  
See instructions for use



Manufactured by:

Péters Surgical  
42, rue Benoît Frachon  
93013 Bobigny Cedex - FRANCE

Tel: +33 (0)1 48 10 62 62  
Fax: +33 (0)1 48 10 22 99  
commercial@peters-surgical.com  
www.peters-surgical.com



REF: IC5 330