


STANDARD GOWNS – TOP

WITH 2 TOWELS AND COVER DRAPE

NON STERILE COMPONENT CODE	STERILE COMPONENT CODE	SIZE	CARTON
S04ALL0052TOP	F21CAM0052TOP	M	40 (20X2)
S04ALL0053TOP	F21CAM0053TOP	L	36 (18X2)
S04ALL0054TOP	F21CAM0054TOP	XL	36 (18X2)
S04ALL0056TOP	F21CAM0056TOP	XXL	36 (18X2)



SIZE	BODY DIMENSIONS		ARM DIMENSIONS	
	LENGHT (cm)	BREADTH (cm)	LENGHT (cm)	LENGHT (cm)
M	120	132	58,5	57
L	135	141	59,5	57
XL	145	155	61	63
XXL	155	155	61	70

FEATURES

• ALCOHOL RESISTANT	• ANTI-INFECTION CARD
• VELCRO NECK CLOSURE	• REINFORCED AND ULTRASONIC BONDED BELTS
• YELLOW SMMS EDGES	• 2 TOWELS
• ARMS ULTRASONIC BONDED	• DOUBLE WRAPPED
• CUFFS IN TUBULAR JERSEY	• 2 TRACEABILITY LABELS

MATERIALS: Unimat 1 – 13 – 22 – 23

PACKAGING: PAPER PE POUCH

CARTON SIZE: 60X40X38 cm

STERILIZATION: E/O VALIDITY 5 YEARS

STORAGE: KEEP IN A COOL AND DRY PLACE

CEE 93/42 DIRECTIVE CLASSIFICATION: CLASS I STERILE

CND CLASSIFICATION: T020401

REPERTORY: 456047

Latex Free